

APPENDIX A

CLEANING SAMPLES FOR CD/CA

- 1) Prepare Samples
 - a. Label 0.65 mL centrifuge tubes with sample ID—top and side
 - b. Weigh samples ~0.5+mg
 - i. Tare centrifuge tube and add samples with tweezers (stainless)
- 2) Precleaning
 - a. Water Wash
 - i. add ~250 μL mQ-water
 - ii. ultrasonicate 20 minutes
 - iii. resuspend
 - iv. repeat until solution is clear
 - b. Oxidizing Solution
 - i. add 1 x 250 μL 50:50 mixture of 30% H_2O_2 /2M NaOH
 - ii. ultrasonicate 20 min
 - iii. siphon off solution
 - iv. resuspend 3 times with mQ-water
 - c. Water Rinse
 - i. add ~250 μL mQ-water
 - ii. ultrasonicate 2 minutes
 - iii. resuspend
 - d. Methanol
 - i. add 1 x 250 μL CH_3OH
 - ii. ultrasonicate 5 minutes
 - iii. resuspend 3 times with mQ water
 - e. Water Rinse
 - f. Bleach
 - i. add 1 x 250 μL 50:50 mixture of 30% H_2O_2 /1% HClO_4
 - ii. wait 30–60 s
 - iii. resuspend 2–3 times with mQ-water
 - g. Water Rinse
 - h. Transfer to acid leached polyethylene vials for cleaning
- 3) Cleaning
 - a. Water Rinse—3 times
 - b. Methanol Wash
 - i. add 1 x 250 μL CH_3OH
 - ii. ultrasonicate 2 min
 - iii. resuspend 2–3 times with mQ-water
 - c. Water Rinse

- d. Nitric Acid Wash
 - i. add 100 μL 0.2% HNO_3
 - ii. ultrasonicate 1 minute
 - iii. quickly check pH (use 10 μL pipet)
 - iv. resuspend 2–3 times with mQ-water
 - e. Water Rinse
 - f. Oxidizing Solution
 - i. add 1 x 250 μL 50:50 mixture of 30% H_2O_2 /0.1M NaOH
 - ii. cover rack tightly
 - iii. heat for 20 min in boiling water bath, ultrasonicate 5–10 s every 2 min (starting at 0), invert after each ultrasonication
 - iv. Rinse 2–3 times with mQ-water
 - g. Water Rinse–2 times
 - h. Reducing Solution–ALL IN HOOD
 - i. Make up working solution
 - 1. ~15 mL NH_3 /citrate solution
 - 2. ~15 mL fresh NH_4OH
 - 3. 1 x 1000 μL anhydrous hydrazine
 - ii. Add 1 x 250 μL working solution
 - iii. Heat for 30 min in boiling water bath, ultrasonicate 5–10 s every 2 min (starting at 0), invert after each ultrasonication
 - iv. Resuspend 2–3 times with mQ-water
 - i. Water Rinse–2 times
 - j. Oxidizing Solution
 - k. Water Rinse–2 times
 - l. Nitric Acid Wash
 - m. Water Rinse–2 times
- 4) Transfer and Dissolve
- a. Transfer
 - i. Add 1 x 50 μL mQ-water to each tube
 - ii. Use 10 μL pipet with cut off pipet tip to transfer sample with a minimum of solution to a new labeled tube
 - b. Postcleaning
 - i. Add 1 x 100 μL 0.2% HNO_3
 - ii. Leave for ~1 min
 - iii. Aspirate
 - iv. Resuspend 2–3 times with mQ-water
 - c. Dissolve
 - i. Add 200/250 μL 2% HNO_3 and wait for sample to dissolve